
PHE National Parasitology Reference Laboratory, Hospital for Tropical Diseases, 3rd Floor Mortimer Market, Centre, Capper Street, London WC1E 6JB, TEL: +44 (0) 207 383 0482, FAX +44 (0) 207 388 8985

Modified Ziehl-Neelsen


Use of the modified Ziehl-Neelsen stain for faecal smears has already been established for coccidian protozoa, in particular, oocysts of *Cryptosporidium* species, but it is also useful to confirm the presence of oocysts of *Isospora belli* and *Cyclospora cayetanensis*. (See Colour Plate 3, page 55)

Method

- a) Faecal smears are made either directly from the stool sample or from the concentration deposit.
- b) Allow to air dry.
- c) Fix in methanol for 3 minutes.
- d) Stain with strong carbol fuchsin for 15-20 minutes.
- e) Rinse thoroughly in tap water.
- f) Decolourise in acid alcohol (1% HCl in methanol) for 15-20 seconds.
- g) Rinse thoroughly in tap water.
- h) Counterstain with 0.4% malachite green (or methylene blue) for 30-60 seconds.
- i) Rinse thoroughly and air dry.
- j) Examine using x40 and x100 objectives.


Oocysts of *Cryptosporidium parvum*

Oocysts of *Cryptosporidium parvum* do not concentrate well using standard concentration techniques and are identified using various staining techniques. Using the Modified Ziehl-Neelsen stain, the oocysts are acid-fast. However, staining within a smear and between specimens can vary from unstained, to partial red staining, to complete staining.




Oocysts of *Isospora belli*

The oocysts of *Isospora belli* can be demonstrated in faeces after formol-ether concentration. Alternatively, they can be seen in a faecal smear stained by modified Ziehl-Neelsen where the oocysts stain a granular red colour against a green background.



Oocysts of *Cyclospora cayetanensis*

The oocysts of *Cyclospora cayetanensis* can be seen in formol-ether concentrated stool samples. Alternatively they can be seen when stained with modified Ziehl-Neelsen where they exhibit variable staining; some cysts being acid fast whereas others appear as a round hole against the background. Some are seen as glassy wrinkled spheres.



Phenol- auramine stain³

This stain can be used as an alternative to the modified Ziehl-Neelsen stain for staining oocysts of *Cryptosporidium parvum*.

Method

- a) Make faecal smears as for ZN and fix in methanol.
- b) Stain with phenol-auramine (Lemperts) for 10-15 minutes.
- c) Rinse thoroughly in tap water.

- d) Decolourise in acid alcohol. (as for ZN)
 - e) Rinse thoroughly in tap water.
 - f) Counterstain with 0.1% potassium permanganate for 30 sec.
 - g) Rinse thoroughly in tap water, allow to air dry.
- Do not blot dry, many brands of blotting paper will fluoresce!
- h) Observe the films with blue light under an incident light fluorescent microscope and low power followed by oil immersion x100 objective if oocysts are suspected.

The oocysts of *Cryptosporidium parvum* appear as bright yellow discs against a dark background. The oocysts of *Isospora belli* and *Cyclospora cayetanensis* do not fluoresce well using the phenol-auramine stain.

Summary of Diagnostic Characteristics

Microscopic Characteristics	<i>Cryptosporidium parvum</i>	<i>Isospora belli</i>	<i>Cyclospora cayetanensis</i>
Size	4 - 6 μ	20-33 μ - 10 - 19 μ	8 - 10 μ
Identified in formol-ether concentrate by light microscopy	No	Yes	Yes
Identified by modified Ziehl-Neelsen	Yes	Yes	Yes
Identified by Phenol Auramine stain	Yes	Variable	No